Instruction manual

VWR® Zippette Pro

Zippette Pro 0-30 ml manual priming
Zippette Pro 0-30 ml spring-loaded

North American Catalogue Number: 75856-470

European Catalogue Number: 613-5281


Table of contents

EN Instruction Manual 1 - 14	DE Bedienungsanleitung 29 - 42	IT Manuale d'istruzioni 57 - 70
FR Manuel d'instructions 15 - 28	ES Manual de instrucciones 43 - 56	PT Manual de instruções 71 - 84

United States

Manufactured for:

VWR International, LLC
100 Matsonford rd
Radnor, PA 19087
800 932 5000
us.vwr.com

Europe

Address of legal manufacturer:

VWR Internation bvba
Researchpark Haasrode 2020
Geldenaaksebaan 464
B-3001 Leuven
+ 32 16 385011
be.vwr.com

Contents

WARNING	SPECIFICATION	FLUID-PATH FOR CLEANING OR REPLACEMENT
3	5	10
SAFETY INFORMATION	INSTRUCTIONS FOR USE	STERILISATION/AUTOCLAVING
4	6	11
PACKAGE CONTENTS	PRIMING RESERVOIR	ACCESSORIES & SPARES
4	6	14
INSTALLATION	THE VOLUME ADJUSTMENT KNOB	ADAPTORS
4	6	14
RESTRICTIONS OF USE	USER CALIBRATION	RESERVOIRS
4	6	14
BEFORE USING THE ZIPPETTE PRO	TROUBLESHOOTING	SPARE PARTS
4	6	15
ASSEMBLY	GENERAL MAINTENANCE	TECHNICAL SERVICE
4	6	15
INSERTING	MAINTENANCE/CLEANING	WARRANTY
4	9	15

Warning!

- Always observe general safety regulations (e.g. wear suitable eye, face, hand and general protection).
- Follow this operating instruction manual.
- NEVER use force! Breakage of any part may lead to the operator and/or other persons being exposed to hazardous substances.
- Do not use the instrument where flammable fumes may exist.
- NEVER use damaged or deformed components.

ZIPPETTE PRO

Safety Information

- When using the instrument, take into consideration the chemical compatibility of the materials in contact with the liquid
- Never leave the Zippette Pro on your work bench with the piston barrel full of liquid without securing the dispense anti-drip tap into 'Closed' position
- Always check that the dispense anti-drip tap is in the 'Open' position before attempting to dispense
- Before adjusting volume always place receiving vessel under nozzle and open tap
- Never use force
- Use only original spare parts and original accessories (See pages 14-15).
- If in doubt, consult your safety officer

Package contents

The Zippette Pro kit comprises:

- Zippette Pro dispenser
- Adaptors: 38 mm, 40 mm and 45 mm
- PTFE inlet tube
- Return spring

Installation

Restrictions of Use

NEVER use the Zippette Pro with:

- Liquids which are not compatible with PTFE, PVDF, FEP, borosilicate glass, alumina ceramic or those which might attack platinum-iridium
- Hydrofluoric acid
- Liquids which contain solid particles
- Temperature limits of Zippette Pro and reagent are +15 °C to +40 °C

Before using the Zippette Pro



Check that the instrument has not been damaged in transit.

Assembly

The Zippette Pro is packed with the dispense nozzle attached and the inlet feed tube removed. The 300 mm length of PTFE inlet tube provided should be trimmed to fit your particular reservoir. The bottom end should be cut diagonally at a length that is close to the bottom of the container when the Zippette Pro is assembled to the container. Longer lengths of inlet tube are available on request.

The threaded platform base of the Zippette Pro has a 30 mm screw thread for use with the recommended accessories displayed on page 14. The assembled Zippette Pro is screwed to the reservoir using gentle hand torque applied to the Zippette Pro's threaded platform base only. Removal should also be by means of hand torque applied to this same base.

Do not operate the piston until the unit is safely and fully mounted on the reservoir bottle.


Inserting the return spring (when applicable)


To insert the return self-priming spring follow the steps below:

1. Make sure that the Zippette Pro is completely empty. For ease of operation, assemble the Zippette Pro onto an empty reservoir bottle and turn the dispense anti-drip tap to 'Open' position.
2. Turn the volume adjustment knob, to set the Zippette Pro to the max dispense volume of 30 ml (Drawing 1).
3. Turn the fixing ring which is located at the bottom of the Zippette Pro's piston sleeve until its raised rib is aligned with the front piston sleeve slot and remove by gently pulling down (Drawing 1). The fixing ring should stay on the barrel protection inner sleeve (Drawing 2).
4. Press the two fixing levers on the back of the piston sleeve inwards (Drawing 2), and lever the back piston sleeve away from the front sleeve assembly (Drawing 3).
5. Remove the front sleeve assembly from the piston locator and move the piston out of its glass barrel (Drawing 3 & 4).
6. Carefully guide the piston through the spring, and locate the piston back into its glass barrel (Drawing 5).
7. Reassemble the front Zippette Pro piston sleeve (Drawing 5) making sure that the cursor locates under the stop ring (Drawing 5).
8. Refit the back piston sleeve and pull up the fixing ring (Drawing 6) – turn it half circle to locate it onto the Zippette Pro piston sleeve (Drawing 7). The unit is now ready for use (see operating instructions page 7).

Drawing 1


Drawing 2


ZIPPETTE PRO


Drawing 3


Drawing 4


Drawing 5

Drawing 6

Drawing 7

Specification

Accuracy $\leq \pm 0.3\%$ on maximum delivery and a precision $\leq 0.1\%CV$ using distilled water at 20 °C.

Testing takes place according to DIN EN ISO 8655 using maximum nominal volume of distilled water at 20 °C, with a uniform and smooth dispense action.

Instructions for use

Priming reservoir

- Place a container under the Zippette Pro's dispense nozzle.
- Turn the volume adjustment knob to 30 ml to enable free movement of the piston.
- Open the dispense anti-drip tap.
- Prime the unit with a few gentle up and down strokes, taking the piston right down to its lowest stop position and lifting it up approximately 25 mm.
- Repeat until a steady, bubble free flow is visible in the dispense tube.

The volume adjustment knob

Turning the volume adjustment knob (Drawing 1) clockwise, will increase the volume to be dispensed. Each full revolution of the knob represents 2 ml and is shown on the cursor aligned to the graduations.

For fine volume adjustments, turn the knob a quarter turn in either direction for 0.5 ml adjustment.

User calibration procedure

Your Zippette Pro has been laboratory calibrated to its nominal volume of 30 ml. However, due to changes in environmental conditions and the viscosity of the media which you dispense, a recalibration might be required. You can either recalibrate at regular intervals such as once a week or whenever you notice that the dispensed volume is different from the volume displayed by the unit.

To fully recalibrate your Zippette Pro follow the procedure below:

1. Set the Zippette Pro to the nominal volume of 30 ml or any other volume which is the most common volume you dispense.
2. Follow the common rules for calibration used in statistical quality control (ISO 8655-3:2002). Dispense five full volumes of distilled water at 20 °C to establish the actual mean volume of liquid dispensed. If the average result varies from the volume displayed, you should recalibrate the Zippette Pro.
3. Using a Cross Head screw driver no.00, slightly undo the cursor's screw and move the cursor up or down until it aligns with the graduation of the actual volume established in step 2.

Problem	Possible	Remedy
Air bubbles appear in dispense tube	Liquid reservoir is empty	Refill reservoir and prime unit
	Too fast filling action	Pump liquid out, refill and dispense slower
	Glass barrel is not sealing against FEP 'O' ring	Unscrew the threaded platform ring and make sure that the FEP 'O' ring is properly fitted into its recess, if necessary stretch 'O' ring slightly using fingertips
Barrel does not fill with liquid	Leaking piston	Clean piston. If problem persists, replace piston assembly (Drawing 13)
	Leaking discharge valve	Clean unit by flushing through - if problem persists, replace pedestal complete assembly (Drawings 16 & 17)
	Inlet tube not fitted correctly Inlet valve stuck	Connect inlet tube correctly Free inlet valve by inserting a thin rod into the inlet aperture
Dispensing not possible	Glass barrel is not sealing against the FEP 'O' ring	Unscrew the threaded platform ring and make sure that the FEP 'O' ring is properly fitted into its recess, if necessary stretch 'O' ring slightly using fingertips
	Blocked dispense tube nozzle	Disassemble the dispense tube and flush through with cleaning fluid
	Discharge valve stuck	Clean unit by immersing platform in cleaning fluid - if problem persists, replace pedestal complete assembly (Drawings 16 & 17)
Wrong dispense volume	Instrument not calibrated	Calibrate the unit - See page 5 'User Calibration'
	Leaking valves	Clean platform base - if problem persists, replace pedestal complete assembly (Drawings 16 & 17)
Liquid appears between glass barrel and polypropylene protection sleeve	Glass barrel is not sealing against the FEP 'O' ring	Unscrew the threaded platform ring and make sure that the FEP 'O' ring is properly fitted into its recess if necessary stretch 'O' ring slightly using fingertips
	Damaged FEP 'O' ring	Replace FEP 'O' ring

General maintenance

Maintenance/Cleaning

Note: All maintenance should be carried out wearing suitable eye protection and protective clothing. If in doubt, consult your safety officer.

1. Make sure that the Zippette Pro is completely empty and turn the anti-drip tap to 'Open' position.
2. Place the instrument into an empty sink together with its reservoir. Unscrew the threaded platform base from the reservoir and lift the dispenser's intake tube carefully out of the reservoir, whilst tapping it against the reservoir's aperture to shake off any droplets from the intake tube.
3. Hold the dispense nozzle over the aperture of the reservoir and apply gentle piston strokes in order to return any syringe contents into the reservoir.
4. Flush out remaining glass barrel contents with distilled water or a suitable solvent. This will preserve the smooth action of the piston and free action of the inlet and outlet valves. If the inlet valve does stick and is not freed by flushing, it may be freed by gently inserting a thin rod into the inlet aperture and gently pushing the ball off its seating. Re-check unit for operation.
5. Empty the instrument completely after cleaning, and flush through with distilled water.

Fluid-path disassembly procedure for cleaning or replacement

Maintenance/Cleaning

Disassembly should only be undertaken AFTER the unit has been cleaned, using the standard cleaning procedure described on page 9. For autoclaving follow steps 1-3.

1. Undo the piston sleeve fixing ring by turning it to align with the sleeve slot and sliding it down from the sleeve (Drawing 9 & 10).
2. Press the two fixing levers on the back of the piston sleeve inwards, and lever the back piston sleeve away from the front sleeve assembly (Drawing 10 & 11).
3. Remove the front piston sleeve (Drawing 12).
4. Pull out the piston assembly from the glass barrel (Drawing 13).
5. Undo the threaded platform ring of the base platform (Drawing 14) and turn the stop ring clockwise and pull it off the assembly (Drawing 15).
6. Carefully slide off the fixing ring and the threaded platform ring, together with the clear polypropylene protection sleeve (Drawing 16).
7. Take out the glass barrel from the base platform. Under the glass barrel is an encapsulated FEP 'O' ring which should be handled with extreme care (Drawing 16).
8. To remove the FEP dispense tube (Drawing 17): First turn the anti-drip tap to the closed position then lift off the dispense tube protective cover and pull out the dispense tube from its seating.
9. Reassemble the instrument following the above points in reverse order. Make sure that the FEP 'O' ring is well located into the base platform's recess.

Maintenance/Cleaning

- Do not use hard tools to scrape off residue of reagent from the piston
- Take special care not to damage the piston when inserting it into the glass barrel
- Bevelled inside edge of the glass barrel must be at the top end when fitted
- Before tightening the threaded ring fully onto the threaded platform base, rotate the polypropylene protection sleeve so that the tooth on its base locates properly into the notch on the threaded ring
- Make sure the front end of the dispense tube is clipped into the dispense tube protection cover from the underside
- Do not use force in assembly or disassembly
- Always wear protective clothing and eye protection when disassembling or assembling your Zippette Pro. Protective clothing and eye protection should also be worn when dispensing

Sterilisation/Autoclaving

Only the fluid path components of the Zippette Pro are suitable for chemical sterilisation or autoclaving at 121 °C after the normal cleaning procedure has been carried out.

Chemical sterilisation

The fluid path components of the Zippette Pro can be soaked overnight in a dilute (1:1000) solution of sodium hypochlorite. (See page 10 for disassembly procedure).

For soaking the components individually, follow the fluid path disassembly instructions on the previous page.

Chemical sterilisation


Autoclaving of fluid-path components is permissible at 121 °C, 2 bar for 20 minutes, after the normal cleaning procedure has been carried out.

Autoclaving can be carried out only in a part-assembled condition with the piston inserted into the glass barrel (Drawing 8). Follow disassembly instruction (page 10) 1 to 3 (Drawings 9-12). Steam sterilisation of the piston outside the glass barrel may damage it.


Loosen the threaded platform ring (Drawing 8) by turning it anti clockwise two full revolutions. Place the assembly together with the PTFE inlet tube on a suitably soft surface within the autoclave in order to avoid metal contact.

Cool slowly back to ambient temperature after autoclaving and before assembly.


Drawing 8


Drawing 9


Drawing 10


Drawing 11


Drawing 12


Drawing 13


Drawing 14


Drawing 15

Drawing 16


Drawing 17

Accessories and Spares


Adaptors

A range of four adaptors for fitting your dispenser to the reservoir are available along with an angled funnel to aid filling of reservoirs.


33 mm

EU Cat. No. 612-3980
NA Cat. No. 53409-990


38 mm

EU Cat. No. 612-3978
NA Cat. No. NA 53409-992

40 mm

EU Cat. No. 612-3977
NA Cat. No. 53409-994

30 mm

45 mm

EU Cat. No. 612-3979
NA Cat. No. 53409-996

Spare parts

Spare parts	Drawing	EU Cat. No.	NA Cat. No.
FEP Dispense Tube	Drawing 17	613-5668	10782-792
Return Spring	Drawing 5	613-5671	10782-804
Piston Assembly	Drawing 13	613-5672	10782-806
Pedestal Complete Assembly	Drawing 16 & 17	613-5674	10782-810
Fixing Ring	Drawing 9, 10 & 16	613-5675	10782-814
Stop Ring	Drawing 5 & 15	613-5685	10782-834
Protection Sleeve	Drawing 16	613-5687	10782-838
Borosilicate Glass Barrel	Drawing 16	613-5690	10782-844
Encapsulated 'O' Ring	Drawing 16	613-5693	10782-850
Threaded Platform Ring	Drawing 14 & 16	613-5663	10782-852
PTFE Inlet Tube		613-5658	10782-772

See pages 6, 12 and 13 for reference

Technical service

Return for Repair

In the unlikely event of repair, or where the damage to the unit necessitates return to your distributor or the manufacturer, please decontaminate the unit and complete the required decontamination certificate where appropriate. In addition, please give written details of the fault. The decontamination certificate is available from the product ordering page on www.vwr.com.

Warranty

The manufacturer warrants that this product will be free from defects in material and workmanship for a period of three (3) years from date of delivery. If a defect is present, VWR will, at its option and cost, repair, replace, or refund the purchase price of this product to the customer, provided it is returned during the warranty period. This warranty does not apply if the product has been damaged by accident, abuse, misuse, or misapplication, or from ordinary wear and tear. If the required maintenance and inspection services are not performed according to the manuals and any local regulations, such warranty turns invalid, except to the extent, the defect of the product is not due to such non performance.

Items being returned must be insured by the customer against possible damage or loss. This warranty shall be limited to the aforementioned remedies. IT IS EXPRESSLY AGREED THAT THIS WARRANTY WILL BE IN LIEU OF ALL WARRANTIES OF FITNESS AND IN LIEU OF THE WARRANTY OF MERCHANTABILITY.

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IRELAND / NORTHERN IRELAND

VWR International Ltd / VWR
International (Northern Ireland) Ltd
Orion Business Campus
Northwest Business Park
Ballycoolin
Dublin 15
Tel.: +353 (0) 1 88 22 222
sales.ie@vwr.com

ITALY

VWR International S.r.l.
Via San Giusto 85
20153 Milano (MI)
Tel.: +39 02 3320311
info.it@vwr.com

THE NETHERLANDS

VWR International B.V.
Postbus 8198
1005 AD Amsterdam
Tel.: +31 (0) 20 4808 400
info.nl@vwr.com

NORWAY

VWR International AS
Brynsalleen 4,
0667 Oslo
Tel.: +47 22 90 00 00
info.no@vwr.com

POLAND

VWR International Sp. z o.o.
Limbowa 5
80-175 Gdansk
Tel.: +48 58 32 38 200
info.pl@vwr.com

PORTUGAL

VWR International - Material de
Laboratório, Lda
Centro Empresarial de Alfragide
Rua da Indústria, nº 6
2610-088 Amadora
Tel.: +351 21 3600 770
info.pt@vwr.com

SPAIN

VWR International Eurolab S.L.
C/ Tecnología 5-17
A-7 Llinars Park
08450 - Llinars del Vallès
Barcelona
Tel.: +34 902 222 897
info.es@vwr.com

SWEDEN

VWR International AB
Fagerstagatan 18a
163 94 Stockholm
Tel.: +46 (0) 8 621 34 00
kundservice.se@vwr.com

SWITZERLAND

VWR International GmbH
Lerzenstrasse 16/18
8953 Dietikon
Tel.: +41 (0) 44 745 13 13
info.ch@vwr.com

UK

VWR International Ltd
Customer Service Centre
Hunter Boulevard - Magna Park
Lutterworth
Leicestershire
LE17 4XN
Tel.: +44 (0) 800 22 33 44
uksales@vwr.com

CHINA

VWR (Shanghai) Co., Ltd
Bld.No.1, No.3728 Jinke Rd,
Pudong New District
Shanghai, 201203- China
Tel.: 400 821 8006
Email: info_china@vwr.com

INDIA

VWR Lab Products Private Limited
No.139, BDA Industrial Suburb,
6th Main, Tumkur Road, Peenya Post,
Bangalore, India - 560058
Tel.: +91 80 28078400
vwr_india@vwr.com

MIDDLE EAST & AFRICA

VWR International FZ-LLC
DSP Laboratory Complex
125, Floor 01
Dubai, United Arab Emirates
Tel.: +971 4 5573271
info.mea@vwr.com

SINGAPORE

VWR Singapore Pte Ltd
18 Gul Drive
Singapore 629468
Tel.: +65 6505 0760
sales.sg@vwr.com

GO TO VWR.COM FOR THE LATEST NEWS, SPECIAL OFFERS AND DETAILS FROM YOUR LOCAL VWR SUPPORT TEAM